

COMPTECH TAB meeting 2016-1

Jin Zhang

COMPTECH Officer @ APS, ANL

01/22/2016 Argonne, IL

Brief introduction

- X-ray Thermal Diffuse Scattering (TDS)
 - **X-ray thermal diffuse scattering (TDS)** can be used to measure: **single-crystal** elastic properties of any crystalline materials (include **opaque materials**) using regular diffraction setups. Incorporation of **DAC** in TDS experiment is possible.
 - New **anisotropic elasticity measurement method** for high-pressure mineral physics field.
- Universal **Membrane cap** for most DACs
 - convert screw-driven DAC into membrane-driven DAC – **remote pressure control**
 - Universal cap that could **fit different DACs without losing DAC opening**
- Standard “cheap” resistive **heater** for DAC
 - **Cheap, reusable**, standard heaters (~\$10/pc comparing with >\$200/pc)
- **Website** Modification
 - **TAB** meeting ppts, Notes
 - **Tools**: including instructions for downloads
 - **Techniques and Facilities**: into to all COMPRES facilities at APS; Updates to technical development
- **Facilities**
 - PX² commissioning, membrane tubing system in 13BMC
 - Sector 34, add-in viewing and optical system for high-pressure experiments
- **Data base**
 - IEDA-kickoff workshop


Time Line:




X-ray thermal diffuse scattering (TDS)

X-ray thermal diffuse scattering (**TDS**) is caused by lattice thermal vibrations (phonons).

Shape of the TDS signal tells how anisotropic the sound velocities are:



Project timeline:



Major effect during 2015.10-2016.1:

New Data extraction code: 34IDE (python code) and 13IDC, BMC (IDL program within ATREX)

New inversion algorithm: intensity scaling factor fitting within each individual image

X-ray thermal diffuse scattering (TDS)

- Current Status:

- Sample testing:

- Ambient condition: **Si (100) (111)** and **forsterite (111)**
 - High-pressure: **Si (100) ($\sim 6\text{GPa}$)**; **forsterite (111) ($\sim 40\text{GPa}$)**


- Experimental Setup:

- **Angular scan** approach + Pilatus100K:13BMC **13IDC** (beam size $30\mu\text{m}$)
 - **Energy scan** to find Bragg peaks, then monochromatic beam for TDS measurement (MarCCD) **34ID**

- Software development:

- TDS calculation : micro **BV model** for Si and Macro-single crystal **elasticity model** for all crystals
 - **Data inversion** : new optimized inversion code with **better efficiency** for lower symmetry materials
 - **Data extraction** : graphical output from ATREX for Angular scan data (with Przemek)

data output from IgorLaueGo + python for energy scan data (with Ruqing)





34ID

13IDC and BMC

Data extraction: graphical output from ATREX for Angular scan data (with Przemek)


1. Crystal orientation: 2 peaks (partially solved by using >8 peaks) – but SPEC has limitation: no symmetry constrain
2. Detector X-Y to HKL conversion: originally assuming polynomial fit $H(x,y)$ $K(x,y)$, $L(x,y)$




1. Utilizing the obtained 11 peaks and angles to get Lattice parameters with constrained symmetry (SPEC does not do that – it treats everything as triclinic)

2. UB matrix can be read into the ATREX software for hkl calculation of each pixel.


3. Detector calibration can be made according to the mesh scan result.




Data extraction: data output from IgorLaueGo + python code for energy scan data



Make movie: Energy scan for individual peak in ROI



Plot intensity VS image number:



Data extraction: data output from IgorLaueGo + python code for energy scan data

```
In [1]: import numpy as np
import sys
import os
# Input working directory: change folder accordingly
workdir=r'C:\Users\jzhang\Dropbox (UH Mineral Physics)\HIGP_Mineral_Physics_Team\TDS\Fo009GPa34\1DataOriginal'
os.chdir(workdir)
# make sure xmlutils.py is located in current working directory
import xmdutils as xu
print os.getcwd()
```

C:\Users\jzhang\Dropbox (UH Mineral Physics)\HIGP_Mineral_Physics_Team\TDS\Fo009GPa34\1DataOriginal

```
In [2]: ## Define detector info and geometry, for one run, the numbers does not change. ##
## Use Igor to read the saved geometry file.##
mar165 = xu.Detector(2048,2048,162.361/2048,162.361/2048)
mar165_geo = xu.DetectorGeometry([-1.238,-15.689,153.471],
                                [0.00752321,-0.01143958,3.14097174])
```

Detector calibration file

```
In [3]: ## Define crystal orientation (a*,b*,c*) from Igor indexing result ##
astar = np.array((-8.787488, +0.469981, -9.857230)) * 1.0e9
bstar = np.array((+2.723327, +5.083838, -2.185390)) * 1.0e9
cstar = np.array((+6.338666, -5.946492, -5.934293)) * 1.0e9
aas = np.linalg.norm(astar)
bbs = np.linalg.norm(bstar)
ccs = np.linalg.norm(cstar)
samp_rot_mat = np.vstack((astar/aas,bstar/bbs,cstar/ccs)).T
```

Crystal orientation

```
In [5]: ## get q vec of the defined HKL peak ##
filefolder = r'C:\Users\jzhang\Dropbox (UH Mineral Physics)\HIGP_Mineral_Physics_Team\TDS\TDSdata\2015-2\34IDE\BX
filename = r'Escan_54.h5' # use highest image intensity to identify which image the Bragg peak is exactly located
h5img = xu.XmdH5Data(os.path.join(filefolder,filename))

## Calculate the Bragg peak absolute pixel X and Y
# Input fitted Bragg peak positions from Igor, Igor only read relative position
peakX = h5img.roi_startX + 19.4544
peakY = h5img.roi_startY + 20.0588
print "Bragg peak position:", (" ",peakX, peakY,")"

hklname = '(-131)' #Input ideak HKL name
HKLxeng = 23.5985 # Input X-ray energy (keV)
pk_hkl = np.array([-1.0,3.0,1.0]) #Input integer ideal HKL
pk_cen_xyz = xu.pixels_to_xyzs([peakX],[peakY],mar165,mar165_geo)
pk_cen_q = xu.pxys_to_qs(pk_cen_xyz,HKLxeng)
pk_qc = np.linalg.norm(pk_cen_q)
print 'current hkl:',pk_hkl

# test computing hkl of the Bragg peak, obtained numbers should be similar to theoretical HKL
testhkl = xu.qs_to_ortho_hkls(pk_cen_q[0],np.array((aas,bbs,ccs)),samp_rot_mat)
print 'Real hkl is', testhkl
```

Calculate hkl based on
X-ray energy
And pixel position

```
Bragg peak position: ( 662.4544 1060.0588 )
current hkl: [-1.  3.  1.]
Real hkl is [-1.00149768  3.00925398  0.9993432 ]
```

Double check:
ideal hkl
calculated hkl

Data inversion: optimized inversion with **better efficiency** for lower symmetry materials

**In 2015-2
TAB meeting**


Only consider 1st order approximation – single-phonon scattering

Method 1: single crystal elasticity model

Method 2: Born-Von-Karman model


1 model for all, only based on Cij and hkl's.
Change Cij's for different materials.

Micro force constant model for each pair of neighboring atoms – **1 material 1 model**



New code gives exactly the same answer to Si when treating the intensity scaling factors as local variables VS global variables.



Si result: line profiles



Vertical cut


$I_{\text{residual}} \sim 10^2$
 $I_{\text{measure}} \sim 10^4$
Fitting uncertainty
~a few percent

horizontal cut



Forsterite result

beam time in 2016-1:34IDE



34IDE
0GPa

13DC
9GPa

Future work:

1. Figure out what was the problem with forsterite data:

Collaborate with Dr. B. Winkler (Germany):

First principle calculation – predict TDS + Cij's


Treat his theoretically generated TDS as real data and invert the Cij's.

See if this inverted Cij's = Cij's from first principle calculation.

2. Start writing up the first paper describing the method and the our results on single-crystal Si up to ~6 GPa.

Data base – IEDA kickoff workshop

- Importance of having a mineral physics data base:
 - EOS, elasticity, rheology
 - Phase diagram, element partitioning



The screenshot shows the IEDA (Interdisciplinary Earth Data Alliance) website. The header includes the IEDA logo and navigation links: HOME, ABOUT, CONTRIBUTE DATA, SERVICES, DATA COLLECTIONS, COMMUNITY, NEWS, HELP, and CONTACT US. The main content area features two data collection highlights:

PetDB
PetDB - The Petrologic Database gives you interactive access to a global synthesis of chemical, isotopic, and mineralogical analyses. Learn more:

- [Search PetDB](#)
- [Publications that use PetDB](#)
- [Recent data additions](#)

Geochron
Geochron gives access to Geochronology and Thermochronology data of many types. Data are contributed by users into an interactive system for archiving and publishing data through the EarthChem Library.

- [Search Geochron](#)

Contribute Data

IEDA welcomes and encourages investigators to contribute their data to the IEDA, which can then be reused by a diverse community now and in the future.

Sample-based Data

- Analytical geochemistry datasets (rocks, sediments, minerals, fluids)
- Geochemical or petrological syntheses
- Geochronological datasets
- Sample metadata (IGSN registration service)
- Technical reports (analytical methods, data reduction, etc.)

Sensor-based Data

- Derived Geophysical Data (e.g. grids, maps, mosaics)
- Photos and images
- Shipboard, airborne, and terrestrial data collected in the Southern Ocean
- Seismic Reflection Field Data from the academic research community
- Processed Seismic Data
- Technical reports (data reduction, etc.)

Other Data Types

- Experimental datasets
- Software tools (e.g., macros, code, etc.)
- Highlight images & videos, maps, photos, diagrams & schematics

For assistance, please contact info@iedadata.org




What do they have?

Where do we fit to?

Data System	Description	Find Data	Geophysical Data	Oceanographic Data	Geochemical Data	Images	Sample Info
Antarctic and Southern Ocean Data Portal (ASODS)	Data collected in the Southern Ocean	Search	X	X	Links to	X	X
Academic Seismic Portal (ASP) at LDEO	Seismic Reflection Field Data from the academic research community	Search	X				
Academic Seismic Portal (ASP) at UTIG	Processed Seismic Data						
Deep Lithosphere Dataset	Petrology of mantle xenoliths						
EarthChem Portal	Access to geochemical data in federated databases						
EarthChem Library	Geochemical datasets and data syntheses						
Geochron	Geochronological and thermochronological data						
GeoPRISMS Data Portal	Data from interdisciplinary research along continental margins						
Global Multi-Resolution Topography (GMRT)	Global compilation of seafloor bathymetry and land topography						
MARGINS Data Portal	Data from interdisciplinary research along continental margins						
MediaBank	High-quality images and videos for education and outreach						
PetDB	Petrological database of the ocean floor						
Ridge 2000 Data Portal	Data from interdisciplinary research at mid-ocean ridges, from mantle to microbe						
SedDB	Geochemical data of marine sediments						
System for Earth Sample Registration (SESAR)	Sample Catalog & Registry for the International Geo Sample Number						
VentDB	Hydrothermal spring geochemistry data system						
USAP-DCC	Metadata and data from NSF-funded Antarctic research programs						

Contribute Your Data to Ea... x +

[www.earthchem.org/data/contribute](#) Search

HOME ABOUT CONTRIBUTE DATA RESOURCES ACCESS DATA HELP CONTACT US

earthchem library
publication & preservation of data

earthchem portal
single-point access to geochemical databases

data synthesis
PetDB, NAVDAT, SedDB, and other topical data collections


data compliance
data management plans & data compliance reports

contribute data
submit and publish your data in EarthChem data systems

Quick Links
[EarthChem Portal](#)
[EarthChem Library](#)
[PetDB](#)
[NavDat](#)
[SedDB](#)
[Geochron](#)
[Data Templates](#)
[SESAR](#)
[Library of Experimental Phase Relations](#)

Contribute Your Data to EarthChem

EarthChem encourages investigators to contribute their data to the EarthChem data systems so that these data can be discovered and reused now and in the future.



Contribute data to the **EarthChem Library**.

What Types of Data Can Be Contributed?


EarthChem provides repository services for the following data types:

- compositional geochemistry (elemental, isotopic)
- geochronology
- petrography (modal compositions, sample descriptions)
- geochemical synthesis datasets
- experimental datasets (e.g., kinetics data)
- sample metadata
- technical reports
- software tools (e.g., MATLAB® code, EXCEL® macros, etc.)
- images (photos, schematics, maps, etc.)

feedback


Universal Membrane cap for most DACs

- Advantage of Membrane cap – convert screw-driven DAC into membrane-driven DAC
- Remote precise pressure control
- Currently available membrane caps
 - Too many DACs in different sizes – one cap for each type
 - Limitations to the DAC opening – IXS, Brillouin, single x-stal XRD
- Our target:
 - **One cap** for most DACs in different sizes: different spacers
 - **No lose in DAC opening**: ideal for single crystal studies
- **Project timeline:**
 - **Designed** the first edition of membrane cap 01/2015
 - Received all **body parts** of membrane cap 06/2015
 - Stainless steel **tubing** ordered and installed in 13BMC 07-08/2015
 - Gas **membranes** sent back to machine shop for **revision** 09/2015
 - GE digital pressure **controller** ordered 08/2015 received 10/2015
 - GE pressure controller and tubing installed and **user used the system** in 11/2015
 - Received the modified membrane in 11/2015, **tested the membrane caps offline** in 12/2015



Welcome COMPRES users to use the membrane cap + control system from 2016-1!

Membrane system setup at PX²




Users used membrane control and tubing in BMC in 2015-3 run with their own Membrane Caps

Universal Membrane cap for most DACs

Membrane Cap



Regular PC cell in the
membrane Cap




4-pin Cell 100° opening




Membrane



4-pin Cell adapter




4-pin Cell within Membrane Cap



Tested offline in Dec. 2015
successfully

Standard “cheap” resistive heater for DAC


- Commercial W-Al₂O₃ metal ceramics heater can reach 1000 K within 30 s, very stable and cheap (<\$10/pc).
- Ready for use, reproducible specifications, reusable, and not expensive
- Calibration and modification of the W-Al₂O₃ heater: ([procedures online updated](#))



- Project status:
 - 1st test with DAC: with **single** heater – sample chamber reached highest T 450K
 - 1st supply to COMPRES users: (Dan Shim and Kurt Leinenweber)
 - 2nd supply to COMPRES users: (Zhenxian Liu)
 - Other testers: (Bin Chen, Xiaojing Lai, Yi Hu, Jay Bass)

Preferred access to existing beamlines


Sector 34 ID



- Ideal for high P experiments
- Successful high P TDS experiments during 2015-2 and 2015-3 run

Beamline Contact: Wenjun Liu <wjliu@anl.gov>; Ruqing Xu <ruqingxu@anl.gov>; Jon Tischler <tischler@anl.gov>

Portable viewing and ruby system



All light path Integrated into the tube

- Difficult to "misalign"
 - Compact
 - for sector 34IDE
- And other non-high pressure beamlines

<http://comptech.compres.us/>

Thanks! Questions?